


Mulino Frangizolle ad Aspi da Laboratorio

Laboratory Lump Breaker Mill with Milling Discs

MFA Lab

PRODOTTO NUOVO
NEW PRODUCT



- ▶ Manutenzione minima
- ▶ Affidabile e robusto
- ▶ Certificato CE

Mulino Frangizolle ad Aspi da Laboratorio

È impiegato nella frantumazione di materie prime di tipo argilloso e può essere utilizzato nei laboratori grazie alle sue dimensioni ridotte. Le zolle vengono frantumate e ridotte a dimensioni accettabili per le eventuali lavorazioni successive.

Il gruppo macinante è costituito da due rotori su cui sono installati i dischi macinanti (aspi) e sono movimentati da una doppia motorizzazione indipendente.

La macchina è munita di raschiatori per la pulizia dei dischi.



- ▶ Very low maintenance
- ▶ Reliable and sturdy
- ▶ Certified CE

Laboratory Lump Breaker Mill with Milling Discs

It is used to crush wet and dry clayey raw materials and can be used in laboratories due to its small size. The lumps are crushed and reduced to an acceptable size for any subsequent processing.

The milling unit consists of two rotors fitted with the milling discs. The rotors are moved by a separate dual drive unit.

The machine is equipped with scrapers for cleaning the discs.